20000003971

RC 2_112014

IKA®

IKA® RC 2 basic

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This user guide include only the main operations and functions!

You can find and download **the complete operating instructions** in your language on our website (www.ika.com)!

Source language: German

EN

Contents

	Page
EC Declaration of Conformity	7
Explication of warning symbols	7
Safety instructions	7
Correct use	9
Commissioning	9
Menu structure	10

EC Declaration of Conformity

We declare under our sole responsibility that this product corresponds to the regulations 2006/95/EC, 2004/108/EC, 1999/5/EC and 2011/65/EU and conforms with the standards or standardized documents: EN 61010-1, -2-010, -2-051; EN 61326-1 and DIN 12876-1.

Explication of warning symbols

DANGER

Indicates an imminently hazardous situation, which, if not avoided, will result in death, serious injury.

WARNING

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

CAUTION

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

NOTICE

Indicates practices which, if not avoided, can result in equipment damage.

Safety instructions

General information:

- **Read the operating instructions completely before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where it can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health, safety and accident prevention regulations.
- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Do not operate the device in explosive atmospheres, with hazardous substances.
- Protect the device and accessories from bumping and impacting.
- Check the device and accessories for damage before each use. Do not use damaged components.
- Safe operation is only guaranteed with the accessories described in the "**Accessories**" section.
- The socket for the mains cord must be easily accessible.
- Socket must be earthed (protective ground contact).

- The voltage stated on the type plate must correspond to the mains voltage.
- The device can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- Disconnect the mains plug before attaching or changing any accessories.
- Disconnect the mains plug before cleaning, maintenance and transportation of the unit.
- The device must only be opened by trained specialists, even during repair. The device must be unplugged from the power supply before opening. Live parts inside the device may still be live for some time after unplugging from the power supply.

NOTICE

Coverings or parts that can be removed from the device must be put back on the device again to ensure safe operation, for example to keep foreign objects and liquids, etc. from getting into the device.

- The device may only be used as prescribed and as described in these operating instructions. This includes operation by instructed specialist personnel.
- When using critical or hazardous materials in your processes, **IKA®** recommends to use additional appropriate measures to ensure safety in the experiment. For example, users can implement comprehensive monitoring equipment.
- Process pathogenic material only in closed vessels under a suitable fume hood. Please contact **IKA®** application support if you have any question.

WARNING

If the **OFF** switch is not within reach when device is operating, an **EMERGENCY STOP** switch that can be easily accessed must be installed in the work area.

- A laboratory circulator refrigerates and circulates fluid according to specified parameters. This involves hazards due to low temperatures and general hazards due to the device of electrical energy. The user is largely protected by the device of relevant standards. Further hazard sources may arise due to the type of tempering fluid, e.g. by exceeding or undercutting certain temperature thresholds or by the breakage of the container and reaction with the carrier fluid. It is not possible to consider all eventualities. They remain largely subject to the judgment and responsibility of the operator. For this reason, it may become necessary for user to take other precautionary safety measures.
- Insufficient ventilation may result in the formation of explosive mixtures. Only use the device in well ventilated areas.
- Use suitable hoses for connection.
- Secure hoses and tubes against slippage and avoid kinks.
- Check hoses, tubes and bath at regular intervals for possible material fatigue (cracks/leaks).
- When device is used for external circulation, extra precaution must be taken for cold fluid leakage due to damaged hose.

DANGER

Do not start up the device if:

- It is damaged or leaking
- Cable (not only supply cable) is damaged.
- After a power failure during operation, the device may start automatically (depending on operating mode).
- Transport the device with care.
- Do not transport or empty the bath while it is still cold. Check the temperature of the bath fluid before draining.
- Always empty the bath before moving the unit.

NOTICE

Always empty the bath if the unit is unused for long time.

Fluids:

CAUTION

Only use the fluids, which fulfill the requirements for safety, health and device compatibility.

Be aware of the chemical hazards that may be associated with the bath fluid used. Observe all safety warning for the fluids.

- Depending on the bath fluid used and the type of operation or toxic can arise. Ensure suitable extraction.
- Do not use any fluid which may cause dangerous reactions during processing.
- Only use recommended bath fluid. Only use non-acid and non corroding fluid.

NOTICE

Never operate the device without sufficient fluid! Check the fluid level detection at a regular basis.

- Continuous monitoring of the filling level of the bath fluid is required.
- For optimum temperature stability, the fluids viscosity should be 50 mm²/s or less at its lowest operating temperature. This permits good fluid circulation and minimizes heating from the pump.
- Untreated tap water is not recommended. It is recommended to use distilled water or high purity water (ion exchangers) and add 0.1 g soda (sodium carbonate Na₂CO₃) /liter, to reduce corrosive properties.

NOTICE

Don't use following fluids:

- Untreated tap water
- Acids or bases
- Solutions with halides: chlorides, fluorides, bromides, iodides or sulfur
- Bleach (Sodium Hypochlorite)
- Solution with chromates or chromium salts
- Glycerine
- Ferrous water.

Correct use

• Use:

Use **RC** (Refrigerated **C**irculators) for cooling and circulating fluids.

Intended Use: Tabletop device

• Range of use (indoor use only):

- Laboratories - Schools
- Pharmacies - Universities

This instrument is suitable for use in all areas except:

- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by **IKA®**.
- If the device is operated improperly or in contrary to the **IKA®** specifications.
- If the device or the printed circuit board are modified by third parties.

Commissioning

NOTICE

After moving the instrument, wait for at least one hour for operation!

Check whether the voltage specified on the type plate matches the mains voltage available.


The power socket used must be earthed (protective earthing).

If above conditions are met, the instrument is ready for operation after plugging in the mains plug.

If these conditions are not met, safe operation is not guaranteed and the device could be damaged.

Observe the ambient conditions (temperature, humidity, etc.) listed under **"Technical Data"**.

After switching on the mains switch (2) at the back of the device or pressing the ON/OFF key (A), all LED segments light up during the self test.

①

Display field

②

Software version (ver)

③

Operating mode (mode)

④

Working setting

After that, following two screens appear alternately.

⑤

⑥

Then the device enters standby status and is ready for operation.

Change the temperature setting with **"Temp (+)"** button (C) or **"Temp (-)"** button (D).

Change the pump speed setting with **"Pump (+)"** button (F) or **"Pump (-)"** button (G).

In standby status, activate the refrigerating function by pressing the “Temp” button (B). The pump and fan starts to run at once. The refrigerating function may delay to start at maximum 2 minutes.

In working status, press the “OK/Pump” button (E) to stop the pump function. The refrigerating function and pump stops.

Note: In running status, the pump speed adjusts automatically according to the load. When the actual speed doesn’t reach the setting speed, the actual speed and setting speed will appear on the screen alternately.

In standby status, press the “OK/Pump” button (E) to start the pump function. The refrigerating function will not be activated.

In working status, press the “Temp” button (B) to stop the refrigerating function, the pump keep running.

The controls of the device allow to be locked by pressing the key button (I), so no accidental changes during operation are possible (LED (J) lights up).

By pressing the key button (I) again, the controls are released (LED (J) no longer lights up).

Menu structure

				Default settings
Menu	MODE	A	-----	Activated
		B	-----	-
		C	-----	-
	Maximum temperature (HI T) -----			40 °C
	Minimum temperature (LO T) -----			-20 °C
	Maximum speed (HI R) -----			3200 rpm
	Minimum speed (LO R) -----			1000 rpm
	Fluid type (FLUI) -----			8
	The type of temperature PID control arithmetic (AUTO)	AUTO 1	-----	Activated
		AUTO 0	Proportional coefficient of PID (Kp 1) --	3.0
			Integrate time of PID (Ti 1) -----	20.0
			Differential time of PID (Td 1) -----	2.0
	Alarm and key tone (BEEP)	BEEP 0	-----	-
		BEEP 1	-----	Activated
	Calibration (CALI)	CALI 0	-----	Activated
		CALI 2	-----	-
		CALI 3	-----	-

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IKA®-Werke GmbH & Co.KG